

HCG antibody

Product Information

| | |
|---------------|---|
| Catalog No.: | FNab10818 |
| Size: | 1mg |
| Form: | liquid |
| Purification: | protein A+G purified |
| Purity: | ≥95% as determined by SDS-PAGE |
| Host: | Mouse |
| Clonality: | monoclonal |
| Clone ID: | 5D1 |
| IsoType: | IgG1 |
| Storage: | PBS with 0.1% proclin300, pH 7.4, -20°C for 12 months(Avoid repeated freeze / thaw cycles.) |

Background

None

Immunogen information

| | |
|--------------|--|
| Immunogen: | purified Chorionic gonadotropin |
| Synonyms: | Choriogonadotropin subunit beta 3 Choriogonadotropin subunit beta (CG-beta) Chorionic gonadotropin chain beta CGB3 CGB |
| Observed MW: | None |
| Uniprot ID : | P0DN86 |

Application

| | |
|-----------------------|--|
| Reactivity: | Human |
| Tested Application: | ELISA, CLIA, GICA(capture antibody) |
| Recommended dilution: | Each laboratory should determine an optimum working titer for use in its particular application. |