

TNFRSF10A antibody

Product Information


Catalog No.:	FNab10828
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	tumor necrosis factor receptor superfamily, member 10a
Synonyms:	Tumor necrosis factor receptor superfamily member 10A Death receptor 4 TNF-related apoptosis-inducing ligand receptor 1 (TRAIL receptor 1, TRAIL-R1) TNFRSF10A APO2 DR4 TRAILR1
Observed MW:	44-50 kDa
Uniprot ID :	O00220

Application

Reactivity:	Human, Rat
Tested Application:	ELISA, WB, IF
Recommended dilution:	WB: 1:500-1:2000; IF: 1:20-1:200
Image:	



Rat lung tissue were subjected to SDS PAGE followed by western blot with FNab10828(DR4 antibody) at dilution of 1:1000