

PCT antibody

Product Information

Catalog No.:	FNab10901
Size:	100µg
Form:	liquid
Purification:	protein A+G purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	8G5
IsoType:	IgG1
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

Background

None

Immunogen information

Immunogen:	pro-Calcitonin
Synonyms:	Calcitonin Calcitonin Katacalcin Alternative names: Calcitonin carboxyl-terminal peptide (CCP) PDN-21 CALCA CALC1
Observed MW:	15 kDa
Uniprot ID :	P01258

Application

Reactivity:	Human
Tested Application:	ELISA, Immunofluorescence chromatography