

## GPR78 antibody

### Product Information


Catalog No.:	FNab10915
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	G-protein coupled receptor 78
Synonyms:	G-protein coupled receptor 78 GPR78
Observed MW:	39 kDa
Uniprot ID :	Q96P69

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500-1:2000
Image:	



Rat kidney were subjected to SDS PAGE followed by western blot with FNab10915 at dilution of 1:1000.