
anti- A29L antibody

Product Information

Catalog No.:	FNab10955
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Immunogen:	Monkeypox Virus (MPXV) Protein A29L
Uniprot ID :	Q77HM6

Application

Reactivity:	Monkeypox Virus
Tested Application:	ELISA, WB
Recommended dilution:	WB: 1:1000-1:5000; ELISA: 1:5000 - 1:10000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com