

## ISL1 antibody

### Product Information


Catalog No.:	FNab10975
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
Clone ID:	None
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months(Avoid repeated freeze / thaw cycles.)

### Immunogen information

Immunogen:	Insulin gene enhancer protein ISL-1
Synonyms:	ISL1 Islet-1 ISLET1 ISL LIM homeobox 1 Isl 1
Observed MW:	50 kDa
Uniprot ID :	P61371

### Application

Reactivity:	Human
Tested Application:	WB
Recommended dilution:	WB: 1:500-1:2000
Image:	



HeLa cell lysates were subjected to SDS PAGE followed by western blot with FNab10975 (ISL1 antibody) at dilution of 1:500