

## Phospho-JAK1(Y1022/Y1034) antibody

### Product Information

Catalog No.:	FNab11002
Size:	100µg
Form:	liquid
Purification:	Immunogen affinity purified
Purity:	≥95% as determined by SDS-PAGE
Host:	Rabbit
Clonality:	polyclonal
IsoType:	IgG
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Background


Phospho-JAK1 (Tyr1022/Tyr1034] Antibody detects endogenous levels of JAK1 only when phosphorylated at Tyr1034, which site historically referenced as Tyr1022.

### Immunogen information

Immunogen:	Peptide
Synonyms:	Tyrosine-protein kinase JAK1 Janus kinase 1 (JAK-1) JAK1 JAK1A JAK1B
Observed MW:	130 kDa
Uniprot ID :	P23458

### Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution: WB:	1:500 - 1:2000
Image:	



3T3-L1 cells treated by SS18H+LIF(m) (25 ng/ml) for 30 minutes were subjected to SDS PAGE followed by western blot with FNab11002 at dilution of 1:1000