

## **FineTest®Violet 610 Anti-Mouse Ly6G/Ly6C (Gr-1) Antibody(RB6-8C5)**

### **Product Information**

Catalog No.:	FV610-30143
Conjugation:	FineTest®Violet 610
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	RB6-8C5
IsoType:	IgG2b, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### **Immunogen information**

Synonyms:	Gr-1 Gr1 Ly-6G/Ly-6C Ly6G/Ly6C Ly6G Ly6C
Uniprot ID :	P35461 P0CW03

### **Application**

Reactivity:	Mouse
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image: