

HRP-conjugated Beta Actin Monoclonal antibody

Product Information

Catalog No.:	HRP-FNab00873
Size:	100µl
Form:	liquid
Purification:	Protein A+G purification
Purity:	≥95% as determined by SDS-PAGE
Host:	Mouse
Clonality:	monoclonal
Clone ID:	0A0
IsoType:	IgG2a
Storage:	PBS with 0.1% proclin300, 50% glycerol and 1%BSA, pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Background

This gene encodes one of six different actin proteins. Actins are highly conserved proteins that are involved in cell motility, structure, and integrity. This actin is a major constituent of the contractile apparatus and one of the two nonmuscle cytoskeletal actins.

Immunogen information

Immunogen:	actin, beta
Synonyms:	ACTB, Actin, actin, beta, Actin, cytoplasmic 1, B actin, Beta actin, F actin, β actin
Observed MW:	42 kDa
Uniprot ID :	P60709

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	ELISA, WB
Recommended dilution:WB:	1:1000 - 1:10000

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com