
PerCP-Cy5.5 Mouse IgG2b, κ Isotype Control(MPC-11)

Product Information

Catalog No.:	PC55-30078
Conjugation:	PerCP-Cy5.5
Size:	20T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	MPC-11
IsoType:	IgG2b, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Synonyms:	IgG2b
Uniprot ID :	None

Application

Reactivity:	None
Tested Application:	FC

Recommended dilution: Volume per test: 5 μ L. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 μ L of antibody per test (million cells in 100 μ L staining volume or per 100 μ L of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image: