

PE Anti-Mouse CD117 (2B8)

Product Information

Catalog No.:	PE-30136
Conjugation:	PE
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	2B8
IsoType:	IgG2b, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Synonyms:	CD117 Kit Mast/stem cell growth factor receptor Kit Proto-oncogene c-Kit SCFR Tyrosine-protein kinase Kit
Uniprot ID :	P05532

Application

Reactivity:	Mouse
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image: