
FineTest®V450 Rat IgG2b, κ Isotype Control(LTF-2)

Product Information

| | |
|--------------|--|
| Catalog No.: | VIO-30103 |
| Conjugation: | FineTest®Violet 450 |
| Size: | 20T/100T |
| Form: | liquid |
| Host: | Rat |
| Clonality: | monoclonal |
| Clone ID: | LTF-2 |
| IsoType: | IgG2b, kappa |
| Storage: | PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Immunogen information

| | |
|--------------|-------|
| Synonyms: | IgG2b |
| Uniprot ID : | None |

Application

| | |
|---------------------|------|
| Reactivity: | None |
| Tested Application: | FC |

Recommended dilution: Volume per test: 5 μ L. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 μ L of antibody per test (million cells in 100 μ L staining volume or per 100 μ L of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image: