

HRP-conjugated Streptavidin

Catalogue No.:	FNSA-0001
Size	100µl
Conjugation:	HRP purified from horseradish root
Formulation:	0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, pH7.6
Application:	ELISA, WB
Suggested dilution:	1:2000-1:10000 for western blotting with ECL substrates 1:1000-1:20000 for ELISA and WB with chromogenic substrates
Storage:	2-8°C for 12 months
Safety Note:	This product is intended for research and manufacturing uses only. It is not a diagnostic device. Product degradation will result from multiple freeze/thaw cycles. It is suggested that the antigen be stored in use size aliquots and thawed just prior to use. This material has been inactivated, however as with all biological materials, it should be handled as potentially infectious. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.