

种属	货号	产品名称	灵敏度	检测范围	目录价 96T
Human	AQ-H0201	Human IL-6 ELISA Kit	0.28pg/ml	0.46-30pg/ml	¥2,880
Bovine	EB0007	Bovine IL-8 ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥3,280
Bovine	EB0016	Bovine PROG ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,480
Bovine	EB0028	Bovine IL-6 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥3,280
Bovine	EB0049	Bovine Testosterone ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,880
Bovine	EB0062	Bovine COR ELISA Kit	0.234ng/ml	0.391-25ng/ml	¥2,480
Bovine	EB0081	Bovine E2 ELISA Kit	7.5pg/ml	12.5-800pg/ml	¥2,480
Bovine	EB0086	Bovine IL-1 β ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥3,280
Bovine	EB0092	Bovine INS ELISA Kit	4.688pg/ml	7.813-500pg/ml(1ng=10uIU)	¥3,280
Bovine	EB0164	Bovine SOD1 ELISA Kit	0.469ng/ml	0.781-50ng/ml	¥3,280
Canine	ECA0002	Canine ADP/Acrp30 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥3,280
Canine	ECA0020	Canine TNF- α ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥3,280
Canine	ECA0042	Canine IL-1 β ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥3,280
Canine	ECA0067	Canine SOD1 ELISA Kit	0.469ng/ml	0.781-50ng/ml	¥3,280
Canine	ECA0103	Canine INS ELISA Kit	4.688pg/ml	7.813-500pg/ml(1ng=10uIU)	¥3,280
Chicken	ECH0061	Chicken PROG ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,480
Chicken	ECH0077	Chicken CORT ELISA Kit	1.688ng/ml	2.813-180ng/ml	¥2,480
Chicken	ECH0113	Chicken SOD1 ELISA Kit	0.469ng/ml	0.78-50ng/ml (20000IU/mg)	¥3,280
Goat	EG0037	Goat IL-1 β ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥3,280
Guinea pig	EGP0005	Guinea pig BDNF ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥3,280
Guinea pig	EGP0020	Guinea pig IgG ELISA Kit	0.938ng/ml	1.563-100ng/ml	¥3,280
Guinea pig	EGP0034	Guinea pig MCP-1 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥3,280
Guinea pig	EGP0049	Guinea pig TNF- α ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥3,280

Guinea pig	EGP0063	Guinea pig ADP/Acrp30 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥3,280
Human	EH0005	Human CXCL1 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Human	EH0008	Human CXCL9 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Human	EH0009	Human EGF ELISA Kit	2.344pg/ml	3.906-250pg/ml	¥2,480
Human	EH0011	Human IL-18 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Human	EH0043	Human BDNF ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Human	EH0073	Human Emmprin ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Human	EH0078	Human CD21 ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,480
Human	EH0090	Human CEA ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,880
Human	EH0099	Human CRP ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Human	EH0110	Human CST3 ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,480
Human	EH0124	Human SELE ELISA Kit	46.875pg/ml	78.125-5000pg/ml	¥2,480
Human	EH0129	Human FGF19 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Human	EH0130	Human FGF21 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Human	EH0151	Human GDNF ELISA Kit	2.813pg/ml	4.688-300pg/ml	¥2,880
Human	EH0152	Human GH ELISA Kit	4.688pg/ml	7.813-500pg/ml	¥2,480
Human	EH0164	Human IFN- γ ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Human	EH0173	Human IL-10 ELISA Kit	4.688pg/ml	7.813-500pg/ml	¥2,480
Human	EH0177	Human IL-15 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Human	EH0184	Human IL-1 α ELISA Kit	2.344pg/ml	3.906-250pg/ml	¥2,480
Human	EH0185	Human IL-1 β ELISA Kit	2.344pg/ml	3.906-250pg/ml	¥2,480
Human	EH0188	Human IL-1RL1 ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,480
Human	EH0199	Human IL-4 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Human	EH0201	Human IL-6 ELISA Kit	2.813pg/ml	4.688-300pg/ml	¥2,480

Human	EH0205	Human IL-8 ELISA Kit	4.688pg/ml	7.813-500pg/ml	¥2,480
Human	EH0210	Human KIM-1 ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,480
Human	EH0214	Human LN ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,880
Human	EH0216	Human LEP ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Human	EH0222	Human MCP-1 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Human	EH0247	Human OPG ELISA Kit	37.5pg/ml	62.5-4000pg/ml	¥2,480
Human	EH0248	Human OPN ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,880
Human	EH0259	Human PRL ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,480
Human	EH0266	Human RBP4 ELISA Kit	0.938ng/ml	1.563-100ng/ml	¥2,880
Human	EH0268	Human REN ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Human	EH0302	Human TNF-α ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Human	EH0311	Human TNFRSF5 ELISA Kit	11.719pg/ml	19.531-1250pg/ml	¥2,480
Human	EH0341	Human PCT ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Human	EH0347	Human MYO ELISA Kit	0.938ng/ml	1.563-100ng/ml	¥2,480
Human	EH0350	Human NT-proBNP ELISA Kit	23.438pg/ml	39.063-2500pg/ml	¥2,480
Human	EH0351	Human CKMB ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,480
Human	EH0359	Human AFP ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,880
Human	EH0363	Human CA19-9 ELISA Kit	1.875IU/ml	3.125-200IU/ml	¥2,480
Human	EH0364	Human CYFRA21-1 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Human	EH0365	Human HE4 ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,880
Human	EH0374	Human INS ELISA Kit	4.688pg/ml	7.813-500pg/ml (1ng=10ul	¥2,480
Human	EH0382	Human FE ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,480
Human	EH0386	Human sTfR1 ELISA Kit	18.75ng/ml	31.25-2000ng/ml	¥2,880
Human	EH0387	Human C-P ELISA Kit	18.75pg/ml	31.2-2000pg/ml	¥2,480

Human	EH0389	Human β -HCG ELISA Kit	4.688mIU/ml	7.813-500mIU/ml (50ng/ml)	¥2,880
Human	EH0393	Human FSH ELISA Kit	0.938mIU/ml	1.563-100mIU/ml	¥2,480
Human	EH0396	Human Ltf ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,880
Human	EH0397	Human LH ELISA Kit	0.563mIU/ml	0.938-60mIU/ml	¥2,480
Human	EH0398	Human CA242 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Human	EH0401	Human TSH ELISA Kit	0.0094uIU/ml	0.0156-1uIU/ml	¥2,480
Human	EH0406	Human MUC1 ELISA Kit	9.375pg/ml	15.625-1000pg/ml(1ng/ml=1000pg/ml)	¥2,480
Human	EH0410	Human GFAP ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,880
Human	EH0415	Human IgA ELISA Kit	1.875ng/ml	3.125-200ng/ml	¥2,480
Human	EH0417	Human IgG ELISA Kit	0.938ng/ml	1.563-100ng/ml	¥2,480
Human	EH0418	Human IgM ELISA Kit	0.938ng/ml	1.563-100ng/ml	¥2,480
Human	EH0420	Human TG ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,880
Human	EH0424	Human B2M ELISA Kit	0.938ng/ml	1.563-100ng/ml	¥2,480
Human	EH0426	Human CA724 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Human	EH0459	Human cTnT ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Human	EH0485	Human IFN- β ELISA Kit	37.5pg/ml	62.5-4000pg/ml	¥2,480
Human	EH0528	Human AMH ELISA Kit	28.125pg/ml	46.875-3000pg/ml	¥2,880
Human	EH0538	Human PAI1 ELISA Kit	0.234ng/ml	0.391-25ng/ml	¥2,480
Human	EH0543	Human S100B ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Human	EH0599	Human SOST ELISA Kit	37.5pg/ml	62.5-4000pg/ml	¥2,880
Human	EH0641	Human COR ELISA Kit	0.234ng/ml	0.391-25ng/ml	¥2,480
Human	EH0649	Human SOD1 ELISA Kit	0.469ng/ml	0.781-50ng/ml(10000IU/mg)	¥2,880
Human	EH0882	Human SHBG ELISA Kit	9.375pmol/L	15.625-1000pmol/L	¥2,880
Human	EH0960	Human PIIINP ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,480

Human	EH1025	Human TTR ELISA Kit	4.688pg/ml	7.813-500pg/ml	¥2,880
Human	EH1092	Human PINP ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,480
Human	EH1396	Human CA50 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Human	EH1397	Human GAST ELISA Kit	0.188pg/ml	0.313-20pg/ml	¥2,880
Human	EH2145	Human DCP ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,880
Human	EH2493	Human AFM ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,880
Human	EH2593-CM	Human ADP/Acrp30 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Human	EH2613	Human ALB ELISA Kit	0.375ng/ml	0.625-40ng/ml	¥2,480
Human	EH2626	Human ANGI ELISA Kit	75pg/ml	125-8000pg/ml	¥2,480
Human	EH2627	Human ANGII ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Human	EH2687	Human BACE1 ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,880
Human	EH2718	Human BNP ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Human	EH2769	Human CC16 ELISA Kit	37.5pg/ml	62.5-4000pg/ml	¥2,880
Human	EH2867	Human COL4 ELISA Kit	0.469ng/ml	0.781-50ng/ml	¥2,480
Human	EH3057	Human FG ELISA Kit	0.938ng/ml	1.563-100ng/ml	¥2,480
Human	EH3082	Human PSA ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,880
Human	EH3267	Human IL-17 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Human	EH3468	Human OC/BGP ELISA Kit	0.141ng/ml	0.234-15ng/ml	¥2,480
Human	EH3654	Human ProGRP ELISA Kit	0.938pg/ml	1.563-100pg/ml	¥2,880
Human	EH3764	Human SELL ELISA Kit	46.875pg/ml	78.125-5000pg/ml	¥2,480
Human	EH4144	Human Ambp ELISA Kit	0.375ng/ml	0.625-40ng/ml	¥2,480
Human	EH4233	Human PGE2 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,880
Human	EH4475	Human PGI ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,880
Human	EH4476	Human PGII ELISA Kit	37.5pg/ml	62.5-4000pg/ml	¥2,880

Human	EH4492	Human TPSA ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,880
Human	EH4850	Human Testosterone ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,880
Human	EH4860	Human CORT ELISA Kit	1.688ng/ml	2.813-180ng/ml	¥2,480
Human	EH9261	Human MT-RNR2 ELISA Kit	0.12pg/ml	0.146-600pg/ml	¥2,880
Hamster	EHA0019	Hamster ADP/Acrp30 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥3,280
Horse	EHS0008	Horse BDNF ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥3,280
Horse	EHS0018	Horse COR ELISA Kit	0.234ng/ml	0.391-25ng/ml	¥2,480
Horse	EHS0025	Horse SOD1 ELISA Kit	0.469ng/ml	0.781-50ng/ml	¥3,280
Mouse	EM0001	Mouse ADP/Acrp30 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Mouse	EM0010	Mouse MPO ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,880
Mouse	EM0020	Mouse BDNF ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Mouse	EM0061	Mouse CRP ELISA Kit	0.234ng/ml	0.391-25ng/ml	¥2,480
Mouse	EM0089	Mouse GM-CSF ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Mouse	EM0094	Mouse IGF1 ELISA Kit	37.5pg/ml	62.5-4000pg/ml	¥2,480
Mouse	EM0109	Mouse IL-1 β ELISA Kit	7.5pg/ml	12.5-800pg/ml	¥2,480
Mouse	EM0121	Mouse IL-6 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Mouse	EM0129	Mouse LEP ELISA Kit	37.5pg/ml	62.5-4000pg/ml	¥2,480
Mouse	EM0135	Mouse MCP-1 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Mouse	EM0160	Mouse SELP ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Mouse	EM0163	Mouse RANK ELISA Kit	37.5pg/ml	62.5-4000pg/ml	¥2,880
Mouse	EM0183	Mouse TNF- α ELISA Kit	2.344pg/ml	3.906-250pg/ml	¥2,480
Mouse	EM0194	Mouse TNFSF11 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Mouse	EM0205	Mouse VEGF ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Mouse	EM0206	Mouse VEGFR3 ELISA Kit	37.5pg/ml	62.5-4000pg/ml	¥2,480

Mouse	EM0210	Mouse IgG ELISA Kit	0.938ng/ml	1.563-100ng/ml	¥2,480
Mouse	EM0254	Mouse FGF21 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Mouse	EM0260	Mouse INS ELISA Kit	46.875pg/ml	78.125-5000pg/ml	¥2,480
Mouse	EM0335	Mouse GFAP ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,880
Mouse	EM0419	Mouse SOD1 ELISA Kit	0.469ng/ml	0.781-50ng/ml(10000IU/mg	¥2,880
Mouse	EM0469	Mouse CC16 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,880
Mouse	EM0481	Mouse PINP ELISA Kit	46.875pg/ml	78.125-5000pg/ml	¥2,480
Mouse	EM0864	Mouse A β 42 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Mouse	EM0946	Mouse CORT ELISA Kit	1.688ng/ml	2.813-180ng/ml	¥2,480
Mouse	EM1155-CM	Mouse IgM ELISA Kit	0.938ng/ml	1.563-100ng/ml	¥2,480
Mouse	EM1204	Mouse MCP-3 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Mouse	EM1245	Mouse NTRK2 ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,880
Mouse	EM1281	Mouse PROG ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,480
Mouse	EM1503	Mouse PGE2 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,880
Mouse	EM1592	Mouse IL-8/cxcl15 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Mouse	EM1618	Mouse TWEAK ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Mouse	EM1721	Mouse COR ELISA Kit	0.234ng/ml	0.391-25ng/ml	¥2,480
Mouse	EM1850	Mouse Testosterone ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,880
Mouse	EM1903	Mouse CD40 ELISA Kit	4.688pg/ml	7.813-500pg/ml	¥2,480
Mouse	EM1905	Mouse LUM ELISA Kit	0.375ng/ml	0.625-40ng/ml	¥2,880
Monkey	EMK0003	Monkey ADP/Acrp30 ELISA Kit	0.938ng/ml	1.563-100ng/ml	¥3,280
Monkey	EMK0060	Monkey PROG ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,480
Monkey	EMK0084	Monkey COR ELISA Kit	0.234ng/ml	0.391-25ng/ml	¥2,480
Monkey	EMK0140	Monkey C-P ELISA Kit	18.75pg/ml	31.2-2000pg/ml	¥3,280

Monkey	EMK0164	Monkey S100B ELISA Kit	4.688pg/ml	7.813-500pg/ml	¥3,280
Monkey	EMK0173	Monkey TTR ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥3,280
Monkey	EMK0189	Monkey TSH ELISA Kit	0.0094uIU/ml	0.0156-1uIU/ml	¥3,280
Porcine	EP0025	Porcine BDNF ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥3,280
Porcine	EP0046	Porcine CST3 ELISA Kit	37.5pg/ml	62.5-4000pg/ml	¥3,280
Porcine	EP0136	Porcine PROG ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,480
Porcine	EP0159	Porcine TNF- α ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥3,280
Porcine	EP0176	Porcine IL-8 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥3,280
Porcine	EP0205	Porcine PGE2 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,880
Porcine	EP0242	Porcine SOD1 ELISA Kit	0.469ng/ml	0.781-50ng/ml	¥3,280
Porcine	EP0254	Porcine COR ELISA Kit	0.234ng/ml	0.391-25ng/ml	¥2,480
Porcine	EP0344	Porcine INS ELISA Kit	4.688pg/ml	7.813-500pg/ml	¥3,280
Rat	ER0003	Rat NGAL ELISA Kit	37.5pg/ml	62.5-4000pg/ml	¥2,880
Rat	ER0006	Rat ADP/Acrp30 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Rat	ER0008	Rat BDNF ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Rat	ER0018	Rat EGF ELISA Kit	2.344pg/ml	3.906-250pg/ml	¥2,480
Rat	ER0020	Rat SELE ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Rat	ER0028	Rat ICAM-1 ELISA Kit	46.875pg/ml	78.125-5000pg/ml	¥2,480
Rat	ER0030	Rat IGF1 ELISA Kit	37.5pg/ml	62.5-4000pg/ml	¥2,480
Rat	ER0042	Rat IL-6 (Interleukin-6) ELISA Kit	37.5pg/ml	62.5-4000pg/ml	¥2,480
Rat	ER0046	Rat LUM ELISA Kit	0.375ng/ml	0.625-40ng/ml	¥2,880
Rat	ER0047	Rat MCP-1 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Rat	ER0048	Rat MCP-3 ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Rat	ER0054	Rat NGF ELISA Kit	23.438pg/ml	39.063-2500pg/ml	¥2,480

Rat	ER0069	Rat VEGF ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Rat	ER0115	Rat LEP ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,480
Rat	ER0229	Rat GFAP ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,880
Rat	ER0332	Rat SOD1 ELISA Kit	0.469ng/ml	0.781-50ng/ml(10000IU/mg	¥2,880
Rat	ER0859	Rat CORT ELISA Kit	1.688ng/ml	2.813-180ng/ml	¥2,480
Rat	ER0950	Rat FGF21 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Rat	ER1014	Rat CXCL1 ELISA Kit	3.75pg/ml	6.25-400pg/ml	¥2,480
Rat	ER1015	Rat Cxcl2 ELISA Kit	4.688pg/ml	7.813-500pg/ml	¥2,480
Rat	ER1048	Rat hs-CRP ELISA Kit	0.094ng/ml	0.156-10ng/ml	¥2,880
Rat	ER1079	Rat IgG ELISA Kit	0.938ng/ml	1.563-100ng/ml	¥2,480
Rat	ER1094	Rat IL-1 β ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Rat	ER1113	Rat INS ELISA Kit	46.875pg/ml	78.125-5000pg/ml	¥2,480
Rat	ER1153	Rat MIP-1 β ELISA Kit	0.938pg/ml	1.563-100pg/ml	¥2,480
Rat	ER1200	Rat NSE ELISA Kit	46.875pg/ml	78.125-5000pg/ml	¥2,480
Rat	ER1255	Rat PROG ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,480
Rat	ER1329	Rat SELL ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Rat	ER1393	Rat TNF- α ELISA Kit	2.344pg/ml	3.906-250pg/ml	¥2,480
Rat	ER1462	Rat Testosterone ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,880
Rat	ER1507	Rat E2 ELISA Kit	7.5pg/ml	12.5-800pg/ml	¥2,480
Rat	ER1626	Rat TWEAK ELISA Kit	9.375pg/ml	15.625-1000pg/ml	¥2,480
Rat	ER1637	Rat ANGII ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Rat	ER1651	Rat COR ELISA Kit	0.234ng/ml	0.391-25ng/ml	¥2,480
Rat	ER1800	Rat PGE2 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,880
Rabbit	ERB0017	Rabbit BDNF ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥3,280

Rabbit	ERB0171	Rabbit IgG ELISA Kit	0.938ng/ml	1.563-100ng/ml	¥3,280
Rabbit	ERB0300	Rabbit PROG ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,480
Rabbit	ERB0301	Rabbit E2 ELISA Kit	7.5pg/ml	12.5-800pg/ml	¥2,480
Sheep	ESH0066	Sheep PGE2 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,880
Sheep	ESH0084	Sheep BDNF ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥3,280
Sheep	ESH0086	Sheep PROG ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,480
Universal	EU0390	Universal E2 ELISA Kit	7.5pg/ml	12.5-800pg/ml	¥2,480
Universal	EU0391	Universal COR ELISA Kit	0.234ng/ml	0.391-25ng/ml	¥2,480
Universal	EU0398	Universal PROG ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,480
Universal	EU0400	Universal Testosterone ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,880
Universal	EU2554	Universal PGE2 ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,880
Universal	EU2580	Universal ALD ELISA Kit	1.875pg/ml	3.125-200pg/ml	¥2,480
Universal	EU2623	Universal ANGII ELISA Kit	18.75pg/ml	31.25-2000pg/ml	¥2,480
Universal	EU3108	Universal CORT ELISA Kit	1.688ng/ml	2.813-180ng/ml	¥2,480
Universal	EU3123	Universal ANGI ELISA Kit	0.188ng/ml	0.313-20ng/ml	¥2,480
Universal	EU3135	Universal OVA ELISA Kit	0.469ng/ml	0.781-50ng/ml	¥2,880