

货号	产品名称	规格
APC-09790	APC Anti-mouse CD150 (TC15-12F12.2)	50T/100T
APC-30001	APC Anti-Mouse CD4 (GK1.5)	50T/100T
APC-30002	APC Anti-Mouse CD3E (145-2C11)	50T/100T
APC-30003	APC Anti-Mouse CD8A (53-6.7)	50T/100T
APC-30004	APC Anti-Human CD3 (OKT-3)	50T/100T
APC-30005	APC Anti-Human CD4 (RPA-T4)	50T/100T
APC-30006	APC Anti-Human CD8A (OKT-8)	50T/100T
APC-30009	APC Anti-Rat CD45 (OX-1)	50T/100T
APC-30010	APC Anti-Rat CD8A (OX-8)	50T/100T
APC-30011	APC Anti-Rat CD3 (G4.18)	50T/100T
APC-30012	APC Anti-Rat CD4 (OX-38)	50T/100T
APC-30013	APC Anti-Mouse CD127(A7R34)	50T/100T
APC-30015	APC Anti-Mouse CD19(1D3)	50T/100T
APC-30016	APC Anti-Mouse CD22(Cy34.1)	50T/100T
APC-30017	APC Anti-Mouse CD25(PC-61.5.3)	50T/100T
APC-30018	APC Anti-Mouse CD279(29F.1A12)	50T/100T
APC-30019	APC Anti-Mouse CD28(37.51)	50T/100T
APC-30020	APC Anti-Mouse CD38(NIMR5)	50T/100T
APC-30021	APC Anti-Mouse CD40(FGK4.5/FGK45)	50T/100T
APC-30022	APC Anti-Mouse CD45(30-F11)	50T/100T
APC-30023	APC Anti-Mouse CD48(HM48-1)	50T/100T
APC-30024	APC Anti-Mouse CD80(16-10A1)	50T/100T
APC-30025	APC Anti-Mouse CD86(GL-1)	50T/100T
APC-30026	APC Anti-Mouse IL-4(11B11)	50T/100T
APC-30027	APC Anti-Mouse Ly6C(Monts1)	50T/100T
APC-30028	APC Anti-Mouse Ly6G(1A8)	50T/100T
APC-30029	APC Anti-Mouse NK1.1(PK136)	50T/100T
APC-30030	APC Anti-Mouse CD274(10F.9G2)	50T/100T
APC-30035	APC Anti-Human CD25(BC96)	50T/100T
APC-30037	APC Anti-Human CD38(HIT2)	50T/100T
APC-30039	APC Anti-Human CD45(HI30)	50T/100T
APC-30048	APC Anti-Human CD11A(TS1/22.1.1.13)	50T/100T
APC-30049	APC Anti-Human CD11B(ICRF44)	50T/100T
APC-30058	APC Anti-Human CD14(M5E2)	50T/100T
APC-30066	APC Anti-Human CD19(CB19)	50T/100T
APC-30072	APC Anti-Mouse/Human CD11b(M1/70)	50T/100T

APC-30073	APC Rat IgG2a, κ Isotype Control(2A3)	20T/100T
APC-30074	APC Anti-Mouse IFN γ (XMG1.2)	50T/100T
APC-30075	APC Rat IgG1, κ Isotype Control(HRPN)	20T/100T
APC-30076	APC Mouse IgG1, κ Isotype Control(MOPC-21)	20T/100T
APC-30077	APC Mouse IgG2a, κ Isotype Control(C1.18.4)	20T/100T
APC-30078	APC Mouse IgG2b, κ Isotype Control(MPC-11)	20T/100T
APC-30083	APC Anti-Human/Mouse CD44 Antibody(IM7)	20T/100T
APC-30084	APC Anti-Mouse CD62L Antibody(Mel14)	50T/100T
APC-30085	APC Anti-Mouse CD11C Antibody(N418)	50T/100T
APC-30087	APC Anti-Rat CD44H Antibody(OX-49)	50T/100T
APC-30092	APC Anti-Human CD45RO Antibody(UCHL1)	20T/100T
APC-30096	APC Anti-Mouse F4/80 Antibody(CI:A3-1)	50T/100T
APC-30097	APC Anti-Mouse CD49b(pan-NK cells) Antibody(DX5)	50T/100T
APC-30098	APC Anti-Mouse IL-17A Antibody(TC11-18H10.1)	50T/100T
APC-30100	APC Anti-Human CD27 (O323)	20T/100T
APC-30103	APC Rat IgG2b, κ Isotype Control(LTF-2)	20T/100T
APC-30105	APC Anti-Mouse CD69 (H1.2F3)	50T/100T
APC-30106	APC Anti-Mouse CD16/32 Antibody(2.4G2)	50T/100T
APC-30109	APC Anti-Human CD69 Antibody(FN50)	20T/100T
APC-30121	APC Anti-Rat CD90/Mouse CD90.1 Antibody[OX-7]	50T/100T
F488-30096	FineTest®488 Anti-Mouse F4/80 Antibody(CI:A3-1)	50T/100T
FITC-09790	FITC Anti-mouse CD150 (TC15-12F12.2)	50T/100T
FITC-30001	FITC Anti-Mouse CD4 (GK1.5)	50T/100T
FITC-30002	FITC Anti-Mouse CD3E (145-2C11)	50T/100T
FITC-30003	FITC Anti-Mouse CD8A (53-6.7)	50T/100T
FITC-30004	FITC Anti-Human CD3 (OKT-3)	50T/100T
FITC-30005	FITC Anti-Human CD4 (RPA-T4)	50T/100T
FITC-30006	FITC Anti-Human CD8A (OKT-8)	50T/100T
FITC-30009	FITC Anti-Rat CD45 (OX-1)	50T/100T
FITC-30010	FITC Anti-Rat CD8A (OX-8)	50T/100T
FITC-30011	FITC Anti-Rat CD3 (G4.18)	50T/100T
FITC-30012	FITC Anti-Rat CD4 (OX-38)	50T/100T
FITC-30013	FITC Anti-Mouse CD127(A7R34)	50T/100T
FITC-30015	FITC Anti-Mouse CD19(1D3)	50T/100T
FITC-30017	FITC Anti-Mouse CD25(PC-61.5.3)	50T/100T
FITC-30018	FITC Anti-Mouse CD279(29F.1A12)	50T/100T
FITC-30020	FITC Anti-Mouse CD38(NIMR5)	50T/100T

FITC-30021	FITC Anti-Mouse CD40(FGK4.5/FGK45)	50T/100T
FITC-30022	FITC Anti-Mouse CD45(30-F11)	50T/100T
FITC-30023	FITC Anti-Mouse CD48(HM48-1)	50T/100T
FITC-30024	FITC Anti-Mouse CD80(16-10A1)	50T/100T
FITC-30025	FITC Anti-Mouse CD86(GL-1)	50T/100T
FITC-30027	FITC Anti-Mouse Ly6C(Monts1)	50T/100T
FITC-30028	FITC Anti-Mouse Ly6G(1A8)	50T/100T
FITC-30029	FITC Anti-Mouse NK1.1(PK136)	50T/100T
FITC-30030	FITC Anti-Mouse CD274(10F.9G2)	50T/100T
FITC-30037	FITC Anti-Human CD38(HIT2)	50T/100T
FITC-30038	FITC Anti-Human CD40(G28.5)	50T/100T
FITC-30039	FITC Anti-Human CD45(HI30)	50T/100T
FITC-30049	FITC Anti-Human CD11B(ICRF44)	50T/100T
FITC-30050	FITC Anti-Human CD11C(BU15)	50T/100T
FITC-30058	FITC Anti-Human CD14(M5E2)	50T/100T
FITC-30066	FITC Anti-Human CD19(CB19)	50T/100T
FITC-30072	FITC Anti-Mouse/Human CD11b(M1/70)	50T/100T
FITC-30073	FITC Rat IgG2a, κ Isotype Control(2A3)	20T/100T
FITC-30074	FITC Anti-Mouse IFN γ (XMG1.2)	50T/100T
FITC-30075	FITC Rat IgG1, κ Isotype Control(HRPN)	20T/100T
FITC-30076	FITC Mouse IgG1, κ Isotype Control(MOPC-21)	20T/100T
FITC-30077	FITC Mouse IgG2a, κ Isotype Control(C1.18.4)	20T/100T
FITC-30078	FITC Mouse IgG2b, κ Isotype Control(MPC-11)	20T/100T
FITC-30084	FITC Anti-Mouse CD62L Antibody(Mel14)	50T/100T
FITC-30085	FITC Anti-Mouse CD11C Antibody(N418)	50T/100T
FITC-30087	FITC Anti-Rat CD44H Antibody(OX-49)	50T/100T
FITC-30092	FITC Anti-Human CD45RO Antibody(UCHL1)	20T/100T
FITC-30096	FITC Anti-Mouse F4/80 Antibody(CI:A3-1)	50T/100T
FITC-30097	FITC Anti-Mouse CD49b(pan-NK cells) Antibody(DX5)	50T/100T
FITC-30100	FITC Anti-Human CD27 (O323)	20T/100T
FITC-30103	FITC Rat IgG2b, κ Isotype Control(LTF-2)	20T/100T
FITC-30105	FITC Anti-Mouse CD69 (H1.2F3)	50T/100T
FITC-30106	FITC Anti-Mouse CD16/32 Antibody(2.4G2)	50T/100T
FITC-30109	FITC Anti-Human CD69 Antibody(FN50)	20T/100T
PE-09790	PE Anti-mouse CD150 (TC15-12F12.2)	50T/100T
PE-30001	PE Anti-Mouse CD4 (GK1.5)	50T/100T
PE-30002	PE Anti-Mouse CD3E (145-2C11)	50T/100T

PE-30003	PE Anti-Mouse CD8A (53-6.7)	50T/100T
PE-30004	PE Anti-Human CD3 (OKT-3)	50T/100T
PE-30005	PE Anti-Human CD4 (RPA-T4)	50T/100T
PE-30006	PE Anti-Human CD8A (OKT-8)	50T/100T
PE-30009	PE Anti-Rat CD45 (OX-1)	50T/100T
PE-30010	PE Anti-Rat CD8A (OX-8)	50T/100T
PE-30011	PE Anti-Rat CD3 (G4.18)	50T/100T
PE-30012	PE Anti-Rat CD4 (OX-38)	50T/100T
PE-30013	PE Anti-Mouse CD127(A7R34)	50T/100T
PE-30015	PE Anti-Mouse CD19(1D3)	50T/100T
PE-30017	PE Anti-Mouse CD25(PC-61.5.3)	50T/100T
PE-30018	PE Anti-Mouse CD279(29F.1A12)	50T/100T
PE-30019	PE Anti-Mouse CD28(37.51)	50T/100T
PE-30020	PE Anti-Mouse CD38(NIMR5)	50T/100T
PE-30021	PE Anti-Mouse CD40(FGK4.5/FGK45)	50T/100T
PE-30022	PE Anti-Mouse CD45(30-F11)	50T/100T
PE-30023	PE Anti-Mouse CD48(HM48-1)	50T/100T
PE-30024	PE Anti-Mouse CD80(16-10A1)	50T/100T
PE-30025	PE Anti-Mouse CD86(GL-1)	50T/100T
PE-30026	PE Anti-Mouse IL-4(11B11)	50T/100T
PE-30027	PE Anti-Mouse Ly6C(Monts1)	50T/100T
PE-30028	PE Anti-Mouse Ly6G(1A8)	50T/100T
PE-30029	PE Anti-Mouse NK1.1(PK136)	50T/100T
PE-30030	PE Anti-Mouse CD274(10F.9G2)	50T/100T
PE-30037	PE Anti-Human CD38(HIT2)	50T/100T
PE-30039	PE Anti-Human CD45(HI30)	50T/100T
PE-30048	PE Anti-Human CD11A(TS1/22.1.1.13)	50T/100T
PE-30049	PE Anti-Human CD11B(ICRF44)	50T/100T
PE-30066	PE Anti-Human CD19(CB19)	50T/100T
PE-30072	PE Anti-Mouse/Human CD11b(M1/70)	50T/100T
PE-30073	PE Rat IgG2a, κ Isotype Control(2A3)	20T/100T
PE-30074	PE Anti-Mouse IFN γ (XMG1.2)	50T/100T
PE-30075	PE Rat IgG1, κ Isotype Control(HRPN)	20T/100T
PE-30076	PE Mouse IgG1, κ Isotype Control(MOPC-21)	20T/100T
PE-30077	PE Mouse IgG2a, κ Isotype Control(C1.18.4)	20T/100T
PE-30078	PE Mouse IgG2b, κ Isotype Control(MPC-11)	20T/100T
PE-30084	PE Anti-Mouse CD62L Antibody(Mel14)	50T/100T

PE-30085	PE Anti-Mouse CD11C Antibody(N418)	50T/100T
PE-30087	PE Anti-Rat CD44H Antibody(OX-49)	50T/100T
PE-30096	PE Anti-Mouse F4/80 Antibody(CI:A3-1)	50T/100T
PE-30097	PE Anti-Mouse CD49b(pan-NK cells) Antibody(DX5)	50T/100T
PE-30103	PE Rat IgG2b,κ Isotype Control(LTF-2)	20T/100T
PE-30104	PE Armenian Hamster IgG Isotype Control[PIP]	20T/100T