

---

## APC Anti-Mouse CD3E(145-2C11)

### Product Information

|              |  |
|--------------|--|
| Catalog No.: | APC-30002  |
| Conjugation: | APC  |
| Size:        | 50T/100T   |
| Form:        | liquid   |
| Host:        | Armenian Hamster   |
| Clonality:   | monoclonal   |
| Clone ID:    | 145-2C11   |
| IsoType:     | IgG  |
| Storage:     | PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months<br>(Avoid repeated freeze / thaw cycles.) |

### Immunogen information

|              |   |
|--------------|---|
| Synonyms:    | T-cell surface glycoprotein CD3 epsilon chain, CD3E, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, CD3E, T3E |
| Uniprot ID : | P22646  |

### Application

|                                       |       |
|---------------------------------------|-------|
| Reactivity:                           | Mouse |
| Tested Application:                   | FC    |
| Recommended dilution:Volume per test: | 5 µ L |