

## APC Anti-Rat CD45 (OX-1)

### Product Information

Catalog No.:	APC-30009
Conjugation:	APC
Size:	50T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	OX-1
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Synonyms:	Receptor-type tyrosine-protein phosphatase C Leukocyte common antigen (L-CA) T200 PTPRC CD45
Uniprot ID :	P04157

### Application

Reactivity:	Rat
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



Rat splenocytes are stained with APC Anti-Rat CD45 Antibody (red histogram). Unstained splenocytes (black histogram) are used as control.