

APC Anti-Mouse CD40(FGK4.5/FGK45)

Product Information

Catalog No.:	APC-30021
Conjugation:	APC
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	FGK4.5/FGK45
IsoType:	IgG2a, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Synonyms:	Tumor necrosis factor receptor superfamily member 5 Cd40 B-cell surface antigen CD40 Bp50 CD40L receptor CD40 Tnfrsf5
Uniprot ID :	P27512

Application

Reactivity:	Mouse
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



KM mouse splenocytes are stained with APC Anti-Mouse CD40 Antibody and PE Anti-Mouse CD19 Antibody.