

APC Anti-Mouse Ly6C(Monts1)

Product Information

Catalog No.:	APC-30027
Conjugation:	APC
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	Monts1
IsoType:	IgG2a, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Synonyms:	Lymphocyte antigen 6 complex locus C Ly 6C Ly6c1 Ly6c2
Uniprot ID :	P0CW03

Application

Reactivity:	Mouse
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



KM mouse bone marrow cells are stained with APC Anti-Mouse Ly6C Antibody (red histogram) or APC Rat IgG2a, κ Isotype Control (black histogram).