

APC Anti-Human CD107A(H4A3)

Product Information

Catalog No.:	APC-30046
Conjugation:	APC
Size:	50T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	H4A3
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Synonyms:	Lysosome-Associated Membrane Protein 1 LGP-120 LAMP-1
Uniprot ID :	P11279

Application

Reactivity:	Human
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



Intracellular staining of the Jurkat cell line with APC Anti-Human CD107a/LAMP-1 Antibody (red histogram) or APC Mouse IgG1 Isotype Control (black histogram).