

## APC Anti-Human IFN $\gamma$ (B27)

### Product Information

Catalog No.:	APC-30053
Conjugation:	APC
Size:	50T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	B27
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Synonyms:	Interferon gamma IFNG IFN-gamma IFN $\gamma$
Uniprot ID :	P01579

### Application

Reactivity:	Human
Tested Application:	FC

Recommended dilution: Volume per test: 5 $\mu$ L. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5  $\mu$ L of antibody per test (million cells in 100  $\mu$ L staining volume or per 100  $\mu$ L of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



HEK293T cells transiently transfected with pcDNA3.1 plasmid encoding Human IFN- $\gamma$  gene are stained with APC Anti-Human IFN- $\gamma$  Antibody (red histogram) or APC Mouse IgG1,  $\kappa$  Isotype Control (black histogram).