

## APC Anti-Human CD9 Antibody(HI9a)

### Product Information

Catalog No.:	APC-30110
Conjugation:	APC
Size:	20T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	HI9a
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Synonyms:	TSPAN29 5H9 antigen CD9 CD9 antigen Cell growth-inhibiting gene 2 protein Leukocyte antigen MIC3 MIC3 MRP-1 Tspan-29 p24
Uniprot ID :	P21926

### Application

Reactivity:	Human
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image: