

APC Anti-Rat CD90/Mouse CD90.1 Antibody(OX-7)

Product Information

Catalog No.:	APC-30126
Conjugation:	APC
Size:	50T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	OX-7
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Synonyms:	Mouse Thy-1.1 Rat Thy-1 Rat CD90 Mouse CD90.1 Thy-1
Uniprot ID :	P01830

Application

Reactivity:	mouse, Rat
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



Rat thymocytes are stained with APC Anti-Rat CD90/Mouse CD90.1 Antibody (red histogram) or APC Mouse IgG1, κ Isotype Control (black histogram).