

## SABC two-component

### Product Information

Catalog No.: ES-0001  
 Size: 200ul/1mL  
 Form: liquid  
 Conjugation: Biotin-HRP and Streptavidin  
 Storage: 2-8°C for 12 months.

### Kit Contents:


Reagents	200 $\mu$ l	1mL
Reagent A	100 $\mu$ L	500 $\mu$ L
Reagent B	100 $\mu$ L	500 $\mu$ L

Preparation of working solutions: Before the tests, Reagent A mix Reagent B on 37°C for 30 Minutes, and use it within 6 hours, the mixture will expire in 12 hours.


### Application

Tested Application: IHC  
 Suggested dilution(mixture): 1:50 for IHC


### Image



Immunohistochemistry of paraffin-embedded mouse lung cancer using PCK antibody and biotin-goat anti rabbit igG, ES-0001(mixture) at dilution of 1:50



Immunohistochemistry of paraffin-embedded Rat Myocardium using  $\alpha$ -SMA antibody and biotin-goat anti rabbit igG, ES-0001(mixture) at dilution of 1:50



Immunohistochemistry of paraffin-embedded mouse breast cancer using PCNA antibody and biotin-goat anti rabbit igG, ES-0001(mixture) at dilution of 1:50

**Note: This product is for research use only, not for diagnostic or therapeutic use.**