

Poly HRP-conjugated Affinipure Goat Anti-Mouse/Rabbit IgG (H+L)


Product Information

CatalogNo.:	ES-0005
Size:	1mL/10mL /100mL
Form:	liquid
Purification:	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Conjugation:	HRP purified from horseradish root
Storage:	PBS with proprietary stabilizer. 2-8°C for 12 months.


Application

Tested Application: IHC, ELISA
Suggested dilution: Ready to use.

Test Results:



Txn (from rabbit IgG) staining in a section of formalin fixed, paraffin embedded human colorectal cancer tissue section



PCNA (from mouse IgG) staining in a section of formalin fixed, paraffin embedded rat liver tissue section

Note: This product is for research use only, not for diagnostic or therapeutic use.