

---

## **Poly HRP conjugated Goat Anti-Rabbit IgG (H+L)--preadsorbed**

### **Product Information**

Catalog No.:	ES-0007
Size:	100 $\mu$ L/500 $\mu$ L/1mL
Form:	liquid
Purification:	The antibody was purified from antisera by immunoaffinity chromatography using rabbit IgG coupled to agarose beads. It is preadsorbed by mouse serum and human serum.
Conjugation:	HRP purified from horseradish root
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months.

### **Application**

Tested Application:	ELISA, WB, IHC
Suggested dilution:	1:5000-1:10000 for western blotting with ECL substrates 1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates 1:50-1:400 for IHC

**Note: This product is for research use only, not for diagnostic or therapeutic use.**