

## **FineTest®488 Anti-Human CD19(CB19)**

### **Product Information**

Catalog No.:	F488-30066
Conjugation:	FineTest®488
Size:	20T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	CB19
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### **Immunogen information**


Synonyms:	B-lymphocyte antigen CD19 Cd19 Differentiation antigen CD19 CD19
Uniprot ID :	P15391

### **Application**

Reactivity:	Human
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



Human peripheral blood lymphocytes are stained with FineTest®488 Anti-Human CD19 Antibody (red histogram). Unstained splenocytes (black histogram) are used as control.