

FineTest®594 Anti-Human VEGFR2/CD309 Antibody

Product Information

Catalog No.:	F594-30163
Conjugation:	FineTest®594
Size:	25T/100T
Form:	liquid
Host:	Rabbit
Clonality:	monoclonal
Clone ID:	3K4
IsoType:	IgG
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Synonyms:	Vascular endothelial growth factor receptor 2 (VEGFR-2) Fetal liver kinase 1 (FLK-1) Kinase insert domain receptor (KDR) Protein-tyrosine kinase receptor flk-1 KDR FLK1 VEGFR2
Uniprot ID :	P35968

Application

Reactivity:	Human
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.