

FineTest®647 Anti-Mouse CD4 (GK1.5)

Product Information


Catalog No.:	F647-30001
Conjugation:	FineTest®647
Size:	25µg/100µg
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	GK1.5
IsoType:	IgG2b, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Synonyms:	T-cell surface glycoprotein CD4 T-cell surface antigen T4 Leu-3 L3T4 T4
Uniprot ID :	P06332

Application

Reactivity:	Mouse
Tested Application:	IF, FC
Recommended dilution:	IF/ICC : 1:250-1:1000; FC:0.1-1 µg/10 ⁶ cells in 100 µL volume.
Image:	



C57BL/6 murine splenocytes are stained with FineTest® 647 Anti-Mouse CD4 Antibody (filled red histogram). Unstained splenocytes (empty black histogram) are used as control.

