

## **FineTest®647 Anti-Mouse TNFA(XT3.11)**

### **Product Information**

Catalog No.:	F647-30166
Conjugation:	FineTest®647
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	XT3.11
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### **Immunogen information**

Synonyms:	Tumor necrosis factor- $\alpha$  Cachectin Necrosin Macrophage cytotoxic factor Differentiation inducing factor TNFSF-2 TNF- $\alpha$  TNF-alpha
Uniprot ID :	P06804

### **Application**

Reactivity:	Mouse
Tested Application:	FC

Recommended dilution: Volume per test: 5 $\mu$ L. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5  $\mu$ L of antibody per test (million cells in 100  $\mu$ L staining volume or per 100  $\mu$ L of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image: