
FineTest®700 Rat IgG2a,κ Isotype Control(2A3)

Product Information

| | |
|--------------|--|
| Catalog No.: | F700-30073 |
| Conjugation: | FineTest®700 |
| Size: | 20T/100T |
| Form: | liquid |
| Host: | Rat |
| Clonality: | monoclonal |
| Clone ID: | 2A3 |
| IsoType: | IgG2a, kappa |
| Storage: | PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.) |

Immunogen information

| | |
|--------------|-------|
| Synonyms: | IgG2a |
| Uniprot ID : | None |

Application

| | |
|---------------------|------|
| Reactivity: | None |
| Tested Application: | FC |

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image: