
FITC Anti-Acetylated tubulin (Lys40)

Product Information

Catalog No.:	FITC-00082
Conjugation:	FITC
Size:	100µL
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	6A10
IsoType:	IgG1
Storage:	PBS with 0.1% proclin300, 50% glycerol and 0.5%BSA, pH 7.3, -20°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Synonyms:	Acetylated alpha tubulin, Acetylated tubulin, acetylated-a-tubulin, acetyl-a-tubulin
Uniprot ID :	P60709

Application

Reactivity:	Human, Mouse, Rat
Tested Application:	IF
Recommended dilution:IF:	1:20-1:200