

---

## FITC Anti-Mouse CD3E (145-2C11)

### Product Information

Catalog No.:	FITC-30002
Conjugation:	FITC
Size:	50T/100T
Form:	liquid
Host:	Armenian Hamster
Clonality:	monoclonal
Clone ID:	145-2C11
IsoType:	IgG
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Synonyms:	T-cell surface glycoprotein CD3 epsilon chain, CD3E, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, CD3E, T3E
Uniprot ID :	P22646

### Application

Reactivity:	Mouse
Tested Application:	FC
Recommended dilution: Volume per test:	5 µ L