
FITC Anti-Mouse CD8A(53-6.7)

Product Information

Catalog No.:	FITC-30003
Conjugation:	FITC
Size:	50T/100T
Form:	liquid
Host:	rat
Clonality:	monoclonal
Clone ID:	53-6.7
IsoType:	IgG2a, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information

Synonyms:	T-cell surface glycoprotein CD8 alpha chain, CD8A, T-lymphocyte differentiation antigen T8/Leu-2, MAL
Uniprot ID :	P01731

Application

Reactivity:	Mouse
Tested Application:	FC
Recommended dilution: Volume per test:	5 µ L