

---

## FITC Anti-Human CD20 (BCA/B20)

### Product Information

Catalog No.:	FITC-30007
Conjugation:	FITC
Size:	50T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	BCA/B20
IsoType:	IgG2a, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information

Synonyms:	B-lymphocyte antigen CD20, MS4A1, B-lymphocyte surface antigen B1, Bp35, Leukocyte surface antigen Leu-16, Membrane-spanning 4-domains subfamily A member 1, MS4A1, CD20
Uniprot ID :	P11836

### Application

Reactivity:	Human
Tested Application:	FC
Recommended dilution: Volume per test:	5 µ L