

## FITC Anti-Rat CD3 (G4.18)

### Product Information

Catalog No.:	FITC-30011
Conjugation:	FITC
Size:	50T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	G4.18
IsoType:	IgG3, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

### Immunogen information


Synonyms:	T-cell surface glycoprotein CD3 $\delta$   $\gamma$   $\epsilon$   and $\zeta$ chains CD3 Complex  T3 CD3
Uniprot ID :	P19377/Q64159/D4A5M2

### Application

Reactivity:	Rat
Tested Application:	FC

Recommended dilution: Volume per test: 5  $\mu$ L. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5  $\mu$ L of antibody per test (million cells in 100  $\mu$ L staining volume or per 100  $\mu$ L of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



Rat splenocytes are stained with FITC Anti-Rat CD3 Antibody (red histogram). Unstained splenocytes (black histogram) are used as control.