

FITC Anti-Human CD127(A019D5)

Product Information

Catalog No.:	FITC-30033
Conjugation:	FITC
Size:	50T/100T
Form:	liquid
Host:	Mouse
Clonality:	monoclonal
Clone ID:	A019D5
IsoType:	IgG1, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Synonyms:	Interleukin-7 receptor subunit alpha IL-7RA CDw127 CD127 IL-7R α
Uniprot ID :	P16871

Application

Reactivity:	Human
Tested Application:	FC

Recommended dilution: Volume per test: 5 μ L. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 μ L of antibody per test (million cells in 100 μ L staining volume or per 100 μ L of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



Human peripheral blood lymphocytes are stained with FITC Anti-Human CD127 Antibody and FineTest®647 Anti-Human CD3 Antibody.