

FITC Anti-Mouse CD14(Sa14-2)

Product Information

Catalog No.:	FITC-30150
Conjugation:	FITC
Size:	50T/100T
Form:	liquid
Host:	Rat
Clonality:	monoclonal
Clone ID:	Sa14-2
IsoType:	IgG2a, kappa
Storage:	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)

Immunogen information


Synonyms:	CD 14 Monocyte differentiation antigen CD14 Myeloid cell-specific leucine-rich glycoprotein
Uniprot ID :	P10810

Application

Reactivity:	Mouse
Tested Application:	FC

Recommended dilution: Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.

Image:



RAW264.7 cells are stained with FITC Anti-Mouse CD14 Antibody (red histogram) or FITC Rat IgG2a Isotype Control (black histogram).