
HRP-conjugated Streptavidin

Product Information

Catalog No.: FNSA-0001
Size: 100µl/500µl
Form: liquid
Conjugation: HRP purified from horseradish root
Storage: 0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300,
pH7.6. 4°C for 24 months.

Application

Tested Application: ELISA, WB

Suggested dilution: 1:2000-1:10000 for western blotting with ECL substrates
1:1000-1:20000 for ELISA and Western blotting with chromogenic
substrates

Note: This product is for research use only, not for diagnostic or therapeutic use.