
HRP-Goat Anti-Rabbit IgG (H+L)

Product Information

Catalog No.: FNSA-0004

Size: 100µl/500µl

Form: liquid

Purification: The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

Conjugation: HRP purified from horseradish root

Storage: 0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6.
-20°C for 24 months.

Application


Tested Application: ELISA, WB, IHC

Suggested dilution: 1:1000-1:10000 for western blotting with ECL substrates


1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates

1:20-1:50 for IHC

Image



Immunohistochemistry of paraffin-embedded human tonsil using CD3E antibody and FNSA-0004(at dilution of 1:20)



Immunohistochemistry of paraffin-embedded mouse brain using VCP antibody and FNSA-0004(at dilution of 1:20)

Note: This product is for research use only, not for diagnostic or therapeutic use.