

---

## HRP-conjugated Affinipure Donkey Anti-Rabbit IgG (H+L)

### Product Information

Catalog No.:	FNSA-0011
Size:	100µl / 500µl
Form:	liquid
Purification:	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Conjugation:	HRP purified from horseradish root
Stabilizer:	3mg/ml BSA ( IgG-Free, Protease-Free)
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

### Application

Tested Application:ELISA, WB

Suggested dilution: 1:1000-1:10000 for western blotting with ECL substrates

1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates

**Note: This product is for research use only, not for diagnostic or therapeutic use.**