

---

## AP-conjugated Affinipure Donkey Anti-Goat IgG (H+L)

### Product Information

|               |   |
|---------------|---|
| Catalog No.:  | FNSA-0022   |
| Size:         | 100µl/500µl   |
| Form:         | liquid  |
| Purification: | The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Conjugation:  | Alkaline Phosphatase purified from calf intestine   |
| Stabilizer:   | 5mg/ml BSA ( IgG-Free, Protease-Free)   |
| Storage:      | 0.01M Tris-HCl, 0.25M NaCl, 50% glycerol, pH8.0.<br>-20°C for 24 months , avoid repeated freezing and thawing.    |

### Application

Tested Application: ELISA, WB

Suggested dilution: 1:1000-1:3000 for ELISA and Western blotting

**Note: This product is for research use only, not for diagnostic or therapeutic use.**