

## **Biotin conjugated Affinipure Donkey Anti-Human IgG (H+L)**

### **Product Information**

Catalog No.:	FNSA-0047
Size:	100µl /500µl
Form:	liquid
Purification:	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Conjugation:	Biotin
Stabilizer:	5mg/ml BSA ( IgG-Free, Protease-Free)
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

### **Application**

Tested Application:	ELISA, WB, IHC/ICC
Suggested dilution:	1:100 – 1:1000 for enzyme immunohistochemistry on tissue sections 1:4000-1:80000 for ELISA and Western blotting with chromogenic substrates

**Note: This product is for research use only, not for diagnostic or therapeutic use.**