
Streptavidin conjugated Affinipure Goat Anti-Human IgG (H+L)

Product Information

Catalog No.:	FNSA-0049
Size:	1mg(1ml)
Form:	liquid
Purification:	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Conjugation:	Streptavidin
Stabilizer:	Antibody Stabilizer (IgG-Free, Protease-Free)
Storage:	4-8°C for 24 months , do not freeze!

Application

Tested Application:	ELISA ,WB
Suggested dilution:	ELISA(1:5000-50000, with biotin-HRP), WB(1:1000-1:10000, with biotin-HRP)

Note: This product is for research use only, not for diagnostic or therapeutic use.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone, Wuhan, Hubei, China(430206)

Tel:(0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com