
Texas Red–conjugated Affinipure Goat Anti-Rabbit IgG (H+L)

Product Information

Catalog No.:	FNSA-0051
Size:	100µl /500µl
Form:	liquid
Purification:	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Fluorophore:	Texas Red sulfonyl chloride (TRSC), Amax=596nm, Emax=620nm
Fluorophore/Protein Absorbance Ratio:	A596/A280=0.57
Stabilizer:	5mg/ml BSA (IgG-Free, Protease-Free)
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

Application

Tested Application: IF, Flow Cyt

Suggested dilution: 1:20 – 1:100 for most applications

Note: This product is for research use only, not for diagnostic or therapeutic use.