
Alexa Fluor 488-Goat Anti-Mouse IgG(H+L)

Product Information

Catalog No.:	FNSA-0055
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified by immunoaffinity chromatography
Conjugation:	Alexa Fluor 488
Fluorophore:	Alexa Fluor 488 carboxylic acid, Amax=493nm; Emax=519nm
Host:	Goat
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

Application

Tested Application:	IF, FC
Recommended dilution:	1:100-1:500

Note: This product is for research use only, not for diagnostic or therapeutic use.