
Rhodamine (TRITC)-Donkey Anti-Goat IgG (H+L)

Product Information

Catalog No.:	FNSA-0067
Size:	100μl/500μl
Form:	liquid
Purification:	The antibody was purified by immunoaffinity chromatography
Fluorophore:	Teramethyl Rhodamine Isothiocyanate(TRITC) Amax=550nm, Emax=570nm
Fluorophore/Protein Absorbance Ratio:	A550/A280=0.49
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

Application

Tested Application:	IF, FC
Recommended dilution:	1:20-1:100

Note: This product is for research use only, not for diagnostic or therapeutic use.