

AMCA conjugated Affinipure Goat Anti-Rabbit IgG (H+L)

Product Information

Catalog No.:	FNSA-0083
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Fluorophore:	7-Amino-4-methylcoumarin-3-ace2c acid A _{max} =350nm, E _{max} =450nm
Fluorophore/Protein Absorbance Ratio:	A ₃₅₀ /A ₂₈₀ =0.76
Stabilizer:	5mg/ml BSA (IgG-Free, Protease-Free)
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6. -20°C for 24 months , avoid repeated freezing and thawing.

Application

Tested Application:	IF, FC
Recommended dilution:	1:20-1:100

Note: This product is for research use only, not for diagnostic or therapeutic use.