

HRP-Goat Anti-Mouse IgM

Product Information

Catalog No.:	FNSA-0091
Size:	100µl/500µl
Form:	liquid
Conjugation:	HRP purified from horseradish root
Purity:	The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
Stabilizer:	3mg/ml BSA (IgG-Free, Protease-Free)
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 0.1% proclin300, 50% glycerol, pH7.6. -20°C for 24 months (Avoid repeated freeze / thaw cycles.).

Application

Tested Application: ELISA, WB

Suggested dilution: 1:2000-1:10000 for western blotting with ECL substrates

1:1000-1:20000 for ELISA and Western blotting with chromogenic substrates

Note: This product is for research use only, not for diagnostic or therapeutic use.