
Alexa Fluor 405-Goat Anti-Mouse IgG(H+L)

Product Information

Catalog No.:	FNSA-0094
Size:	100µl/500µl
Form:	liquid
Purification:	The antibody was purified by immunoaffinity chromatography
Conjugation:	Alexa Fluor 405
Fluorophore:	Alexa Fluor 405 carboxylic acid, Ex: 402nm, Em: 421nm
Host:	Goat
Storage:	0.01M Sodium phosphate, 0.25M NaCl, 30% glycerol, 0.02% Sodium azide, 1%BSA, pH7.6. -20°C for 24 months, avoid repeated freezing and thawing. Store in the dark.

Application

Tested Application:	IF, FC, IHC
Recommended dilution:	1:200-1:2000

Note: This product is for research use only, not for diagnostic or therapeutic use.